

Alaris IV Pump Module

Initial Training Checklist – _____ through _____

These tasks serve as verification of initial hands-on training for the Alaris IV Pump

Name	Unit Date				
Employees Initials	TASK	RN Validator's Initials			
Loading the Alaris Pump module Administration Set					
	Closes IV set roller clamp when safety clamp is open (Prevents free flow				
	when the set is removed from the module)				
Infusion Sets with SmartSite Needle-Free Valves					
	Swabs needle-free valve with a sterile 70% isopropyl alcohol wipe prior				
	to every syringe connection				
	Accesses SmartSite needle-free valve without using blunt cannulas or				
	needles				
	Basic Programming with "Guardrails" Safety Software				
	Articulates what the term "Guardrails" refers to. (Kaiser SCAL drug and				
	fluid library with limits that support best practices.				
	Articulates difference between Hard and Soft Guardrails Limit				
	Selects 'Reprogram'' when Guardrails Hard limit alarms				
	Identifies two (2) types of alerts that are triggered when encounters a hard				
	or soft limit.				
	Demonstrates ability to change value of 'VTBI' after intermittent				
	automatic programming when applicable				
	Demonstrates how to use the Volume/Duration feature				
	Articulates the meaning of a "G" icon showing on the main screen				
	Articulates the meaning of values displayed (Out of Range Limits) when				
	the "G" icon is pressed.				
	Basic Programming <u>without</u> Guardrails Safety Software				
	States when Basic Infusion Programming can be used				
	States difference between 'Basic Infusion' and "Guardrails'				
	programming				
Additional Features					
	Upon infusion completion demonstrates how to use the 'Restore' feature				
	and states what the what the numbers mean				
	Demonstrates how to use the 'Pause' hard key on the Alaris Pump module				

I, the undersigned, verify that I have completed and understand the tasks described above.

Participant's Signature _____

RN Validator's Name (Print)

RN Validator's Signature

Participant successfully programmed pump independently: Y/N _____

Remediation is necessary	because participant was t	unable to perform skills #(s)
(see reverse for list of min		

on this evaluation.

Comments: DLL: C:\Documents and Settings\y538259\Desktop\2012 Final Website Updates 2012\Website Revisions 6.13.2012\Word Docs Format\Alaris Initial Training



Alaris IV Pump Module

Initial Training Checklist – _____ through _

These tasks serve as verification of initial hands-on training for the Alaris IV Pump

Name	Unit	Date
Manager's Name		Notified: Y/N

Minimal qualifying skills competency

Alaris IV Infusion Pump Training

- 1. After repeated tries the participant was not able to correctly program the infusion pump
- 2. Participant bypassed a Soft Guardrail and did not recognize the significance of the action.
- 3. Participant was not able to verbalize what a Soft Guardrail is
- 4. Participant cannot correctly turn the pump on or off

end